



Recombinant DNA-cytosine methyltransferase (dcm)

Product Code	CSB-BP360097EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEE0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQENISVTDS YSTGNAAQAM LEKLLQIYDV KTLVAQLNGV GENHWSAAIL KRALANDSAW HRLSEKEFAH LQTL LKPPA HHPHYAFRFI DLFAGIGGIR RGFESIGGQC VFTSEWNKHA VRTYKANHYC DPATHHFNED IRDITLSHKE GVSDEAAAEH IRQHPEHDV LLAGFPCQPF SLAGVSKKNS LGRAHGACD TQGT LFFDVV RIIDARRPAM FVLENVKNLK SHDQGKTFRI IMQTLDELGY DVADAEDNGP DDPKIIDGKH FLPQHRERIV LVGFRRDLNL KADFTLRDIS ECFPAQRVTL AQLLDPMVEA KYILTPVLWK YLYRYAKKHQ ARGNGFGYGM VYPNNPQSVT RTLSARYYKD GAEILIDRGW DMATGEKDFD DPLNQQHRPR RLTPRECARL MGFEAPGEAK FRIPVSDTQA YRQFGNSVVV PVFAAVAKLL EPKIKQAV AL RQQEAQHGR SR</pre>
Source	Baculovirus
Target Names	dcm
Protein Names	Recommended name: DNA-cytosine methyltransferase EC= 2.1.1.37
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.