



Recombinant Exodeoxyribonuclease 10 (exoX)

Product Code	CSB-EP360128EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEK1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLRIIDTETC GLQGGIVEIA SVDVIDGKIV NPMSHLVRPD RPISPQAMAI HRITEAMVAD KPWIEDVIPH YYGSEWYVAH NASFDRRVLP EMPGEWICTM KLARRLWPGI KYSNMALYKT RKLNVQTPPG LHHHRALYDC YITAALLIDI MNTSGWTAEQ MADITGRPSL MTTFTFGKYR GKAVSDVAER DPGYLRWLFN NLDSMSPELR LTLKHYLENT
Source	E.coli
Target Names	exoX
Protein Names	Recommended name: Exodeoxyribonuclease 10 EC= 3.1.11.- Alternative name(s): Exodeoxyribonuclease X Short name= Exo X Short name= Exonuclease X
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.