



Recombinant Glucarate dehydratase (gudD)

Product Code	CSB-YP360156EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AES3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	SSQFTTPVV TEMQVIPVAG HDSMLMNLGS AHAPFFTRNI VIKDNSGHT GVGEIPGGEK IRKTLEDAIP LVVGKTLGEY KNVLT LVRNT FADR DAGGRG LQTFDLRTTI HVVTGIEAAM LDLLGQHLGV NVASLLGDGQ QRSEVEMLGY LFFVGNRKAT PLPYQSQPDD SCDWYRLRHE EAMTPDAVVR LAEAAEYKYG FNDFKLKGGV LAGEEEAESI VALAQRFPQA RITLDPNGAW SLNEAIKIGK YLKGLSLAYAE DPCGAEQGF S GREVMAEFRR ATGLPTATNM IATDWRQMGH TSLSLQSDIP LADPHFWTMQ GSVRVAQMCH EFGLTWGSHS NNHFDISLAM FTHVAAAAPG KIT AIDTHWI WQEGNQRLTK EPFEIKGGLV QVPEKPLGLV EIDMDQVMKA HELYQKHGLG ARDDAMGMQY LIPGWTFDNK RPCMVR
Source	Yeast
Target Names	gudD
Protein Names	Recommended name: Glucarate dehydratase Short name= GDH Short name= GlucD EC= 4.2.1.40
Expression Region	2-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.