



Recombinant Inosine-guanosine kinase (gsk)

Product Code	CSB-MP360178EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEW8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	MKFPGKRKSK HYFPVNARDP LLQQFQPENE TSAAWVVGID QTLVDIEAKV DDEFIERYGL SAGHSLVIED DVAEALYQEL KQKNLITHQF AGGTIGNTMH NYSVLADDRS VLLGVMCSNI EIGSYAYRYL CNTSSRTDLN YLQGVDPGPIG RCFTLIGESG ERTFAISPGH MNQLRAESIP EDVIAGASAL VLTSYLVRCK PGPEMPEATM KAIEYAKKYN VPVVLTLGTK FVIAENPQWW QQFLKDHVSI LAMNEDEAEA LTGESDPLLA SDKALDWVDL VLCTAGPIGL YMAGFTEDEA KRKTQHPLLP GAIAEFNQYE FSRAMRHKDC QNPLRVYSHI APYMGGPEKI MNTNGAGDGA LAALLHDITA NSYHRSNVPN SSKHKFTWLT YSSLAQVCKY ANRVSQVLN QHSPRLTRGL PEREDSLEES YWDR
Source	Mammalian cell
Target Names	gsk
Protein Names	Recommended name: Inosine-guanosine kinase EC= 2.7.1.73
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.