



Recombinant Escherichia coli Modulator of drug activity B (mdaB)

Product Code	CSB-YP360186ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AEY5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSNILIINGA KKFAHSNGQL NDTLTEVADG TLRDLGHDVR IVRADSDYDV KAEVQNFLWA DVVIWQMPGW WMGAPWTVKK YIDDVFTTEGH GTLYASDGRT RKDPSKKEYGS GGLVQGKKYM LSLTWNAPME AFTEKDQFFH GVGVDGVYLP FHKANQFLGM EPLPTFIAND VIKMPDVPRY TEEYRKHLVE IFG
Source	Yeast
Target Names	mdaB
Protein Names	Recommended name: Modulator of drug activity B
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.