



# Recombinant Human ATP-dependent Clp protease ATP-binding subunit clpX-like, mitochondrial (CLPX)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP005582HU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Uniprot No.</b>       | O76031                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | FTET PAYFASKDGI SKDGSGDGNK KSASEGSSKK SGSGNSGKGG<br>NQLRCPKCGD LCTHVETFVS STRFVKCEKC HHFFVVLSEA DSKKSIIEP<br>ESAAEAVKLA FQQKPPPPPK KIYNYLDKYV VGQSFACKKVL SVAVYNHYKR<br>IYNNIPANLR QQAEVEKQTS LTPRELEIRR REDEYRFTKL LQIAGISPHG<br>NALGASMQQQ VNQQIPQEKR GGEVLDSSHD DIKLEKSNIL LLGPTGSGKT<br>LLAQTLAKCL DVPFAICDCT TLTQAGYVGE DIESVIAKLL QDANYNVEKA<br>QQGIVFLDEV DKIGSVPGIH QLRDVGEGEV QQGLLKLLEG TIVNVPEKNS<br>RKL RGETVQV DTTNILFVAS GAFNGLDRII SRRKNEKYL GFGTPSNLGKG<br>RRAAAAADLA NRSGESNTHQ DIEEKDRLLR HVEARDLIEF GMIPEFVGRL<br>PVVPLHSLD EKTLVQILTE PRNAVIPQYQ ALFSMDKCEL NVTEDALKAI<br>ARLALERKTG ARGLRSIMEK LLEPMFEVP NSDIVCVEVD KEVVEGKKEP<br>GYIRAPTKES SEEEYDSGVE EEGWPRQADA ANS |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Target Names</b>      | CLPX                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: ATP-dependent Clp protease ATP-binding subunit clpX-like, mitochondrial                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Expression Region</b> | 57-633                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                 |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                                                                                                                           |



of lyophilized form is 12 months at -20°C/-80°C.