



Recombinant Human CCN family member 5 (CCN5)

Product Code	CSB-MP026120HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O76076
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	QLCPTPC TCPWPPPRCP LGVPLVLDGC GCCRVCARRL GEPCDQLHVC DASQGLVCQP GAGPGGRGAL CLLAEDDSSC EVNGRLYREG ETFQPHCSIR CRCEDGGFTC VPLCSEDVRL PSWDCPHPRR VEVLGKCCPE WVCGQGGGLG TQPLPAQGPQ FSGLVSSLPP GVPCPEWSTA WGPCSTTCGL GMATRVSNNQN RFCRLETQRR LCLSRPCPPS RGRSPQNSAF
Source	Mammalian cell
Target Names	CCN5
Protein Names	Recommended name: WNT1-inducible-signaling pathway protein 2 Short name= WISP-2 Alternative name(s): CCN family member 5 Connective tissue growth factor-like protein Short name= CTGF-L Connective tissue growth factor-related
Expression Region	24-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.