



Recombinant Nitrogen regulation protein NR (I)

Product Code	CSB-YP360253EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AFB9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MQRGIVVVVD DDSSIRWVLE RALAGAGLTC TTFENGAEVL EALASKTPDV LLSDIRMPGM DGLALLKQIK QRHPMLPVII MTAHSDLDAV VSAYQQGAFD YLPKPFDIDE AVALVERAIS HYQEQQQPRN VQLNGPTTDI IGEAPAMQDV FRIIGRLSRS SISVLINGES GTGKELVAHA LHRHSPRAKA PFIALNMAAI PKDLIESELF GHEKGAFTGA NTIRQGRFEQ ADGGTLFLDE IGDMPLDVQT RLLRVLADGQ FYRVGGYAPV KVDVRIIAAT HQNLEQRVQE GKFREDLFHR LNVIRVHLPP LRERREDIPR LARHFLQVAA RELGVEAKLL HPETEAALTR LAWPGNVRQL ENTCRWLTVM AAGQEVLIQD LPGELFESTV AESTSQMQPD SWATLLAQWA DRALRSGHQV LLSEAQPELE RTLLTTALRH TQGHKQEAAR LLGWGRNTLT RKLKELGME
Source	Yeast
Target Names	glnG
Protein Names	Recommended name: Nitrogen regulation protein NR(I)
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.