



# Recombinant Escherichia coli DNA topoisomerase 4 subunit A (parC)

<b>Product Code</b>	CSB-MP360268ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0AFI2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDMAERLAL HEFTENAYLN YSMYVIMDRA LPFIGDGLKP VQRRIVYAMS ELGLNASAKF KKSARTVGDV LGKYHPHGDS ACYEAMVLMA QPFSYRYPLV DGQGNWGAPD DPKSFAAMRY TESRLSKYSE LLLSELGQGT ADWVPNFDGT LQEPKMLPAR LPNILLNGTT GIAVGMATDI PPHNLREVAQ AAIALIDQPK TTLDQLLDIV QGPDYPTAE IITSRAEIRK IYENGRGSVR MRAVWKKEDG AVVISALPHQ VSGARVLEQI AAQMRNKKLP MVDDLRLDESD HENPTRLVIV PRSNRVDMDQ VMNHLFATTD LEKSYRINLN MIGLDGRPAV KNLLEILSEW LVFRRDTVRR RLNYRLEKVL KRLHILEGLL VAFLNIDEVI EIIIRNEDEPK PALMSRFGLT ETQAEAIEL KLRHLAKLEE MKIRGEQSEL EKERDQLQGI LASERKMNNL LKKEIQADAQ AYGDDRRSPL QEREEAKAMS EHDMLPSEPV TIVLSQMGWV RSAKGHDIDA PGLNYKAGDS FKA AVKGSN QPVVFDSTG RSYAIDPITL PSARGQGEPL TGKLTLP PGA TVDHMLMESD DQKLLMASDA GYGFVCTFND LVARNRAGKA LITLPENAHV MPPVVIEDAS DMLLAITQAG RMLMFPVSDL PQLSKGKGNK IINIPSAEAA RGEDGLAQLY VLPPQSTLTI HVGKRKIKLR PEELQKVTGE RGRRGTLMRG LQRIDRVEID SPRRASSGDS EE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	parC
<b>Protein Names</b>	Recommended name: DNA topoisomerase 4 subunit A EC= 5.99.1.- Alternative name(s): Topoisomerase IV subunit A
<b>Expression Region</b>	1-752
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.