



Recombinant Escherichia coli Spermidine/putrescine-binding periplasmic protein (potD)

Product Code	CSB-MP360281ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AFK9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	DDNNTLY FYNWTEYVPP GLLEQFTKET GIKVIYSTYE SNETMYAKLK TYKDGAYDLV VPSTYYYVDKM RKEGMIQKID KSKLTNFSNL DPDMLNKPFD PNNDYSIPYI WGATAIGVNG DAVDPKSVTS WADLWKPEYK GSLLLTDDAR EVFQMALRKL GYSGNTTDPK EIEAAYNELK KLMPNVAAFN SDNPANPYME GEVNLGMIWN GSAFVARQAG TPIDVWPKE GGIFWMDSLA IPANAKNKEG ALKLINFLLR PDVAKQVAET IGYPTPNLAA RKLLSPEVAN DKTLYPDAET IKNGEWQNDV GAASSIYEEY YQKLGAGR
Source	Mammalian cell
Target Names	potD
Protein Names	Recommended name: Spermidine/putrescine-binding periplasmic protein Short name= SPBP
Expression Region	24-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.