



Recombinant Phage shock protein A (pspA)

Product Code	CSB-BP360291SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AFM8
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	GIFSRFADI VNANINALLE KAEDPQKLVR LMIQEMEDTL VEVIRSTSARA LAEKKQLTRR IEQASAREVE WQEKAELALL KEREDLARAA LIEKQKLTDL IKSLEHEVTL VDDTLARMKK EIGELENKLS ETRARQQALM LRHQAANSSR DVRRLDQSGK LDEAMARFES FERRIDQMEA EAESHSGFKQ KSLDDQFAEL KADDAISEQL AQLKAKMKQD NQ
Source	Baculovirus
Target Names	pspA
Protein Names	Recommended name: Phage shock protein A
Expression Region	2-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.