



Recombinant Bovine Proto-oncogene c-Fos (FOS)

Product Code	CSB-YP008790BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O77628
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MMFSGFNADY EASSSRCSSA SPAGDSLSTYY HSPADSFSSM GSPVNAQDYC TDLAVSSANF IPTVTAISTS PDLQWLQPT LVSSVAPSQT RAPHYGVPT PSAGAYSRAG VMKMTMTGGRA QSIGRRGKVE QLSPEEEEEKR RIRRERNKMA AAKCRNRRRE LTDTLQAETD QLEDEKSALQ TEIANLLKEK EKLEFILAAH RPACKIPDDL GFPEEMSVAS LDLSGGLPEA ATPESSEEAFT LPLLNDPEPK PSVEPVKSVG SMELKAEPFD DYMFPASSRP SGSETARVSP DMDLSGSFYA ADWEPLHGGG LGMGPMATEL EPLCTPVVTC TPSCTTYTSS FVFTYPEADS FPSCAAHRK GSSSNEPSSD SLSSPTLLAL
Source	Yeast
Target Names	FOS
Protein Names	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. In some cases, expression of the FOS gene has also been associated with apoptotic cell death.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.