



Recombinant Corynephage omega Diphtheria toxin, partial

Product Code	CSB-BP360420CQS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00587
Product Type	Recombinant Protein
Immunogen Species	Corynephage omega
Purity	>85% (SDS-PAGE)
Sequence	GADDV VDSSKSFVME NFSSYHGTKP GYVDSIQKGI QKPKSGTQGN YDDDWKGFYS TDNKYDAAGY SVDNENPLSG KAGGVVKVITY PGLTKVLALK VDNAETIKKE LGLSLTEPLM EQVGTEEFIK RFGDGASRVV LSLPFAEGSS SVEYINNWEQ AKALSVELEI NFETRGRGQ DAMYEYMAQA CAGNRVRR
Source	Baculovirus
Protein Names	Recommended name: Diphtheria toxin Short name= DTAAlternative name(s): NAD(+)--diphthamide ADP-ribosyltransferase EC= 2.4.2.36Cleaved into the following 2 chains: 1. Diphtheria toxin fragment A 2. Diphtheria toxin fragment B
Expression Region	26-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.