



# Recombinant Rat Acyl-coenzyme A thioesterase 1 (Acot1)

<b>Product Code</b>	CSB-YP001163RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O88267
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	EATLSLEPA GRSCWDEPLS ITVRGLVPEQ PVTLRAALRD EKGALFRARA LYRADAHGEL DLARAPALGG SFTGLEPMGL IWAMEPERPF WRLVKRDVQT PFVVELEVLD GHEPDGGRLL ARAVHERHFM APGVRRVPVR EGRVRATLFL PPEPGPFGI IDLFGVGGGL LEYRASLLAG KGFAVMALAY YNYDDLPKTM ETMRIEYFEE AVNYLRGHPE VKGPGIGLLG ISKGGELGLA MASFLKGITA AVVINGSVAA VGNTICYKDE TIPPVTILRN QVKMTKDGLK DVVDALQSPL VEQKSFIPVE RSDTTFLFLV GQDDHNWKSE FYANEISKRL QAHGKEKPQI ICYPEAGHYI EPPYFPLCSA GMHLLVGANI TFGGEPKPHS VAQLDAWQQL QTFHFKQLGG KSHGVSPKI
<b>Source</b>	Yeast
<b>Target Names</b>	Acot1
<b>Protein Names</b>	Recommended name: Acyl-coenzyme A thioesterase 1 Short name= Acyl-CoA thioesterase 1 EC= 3.1.2.2 Alternative name(s): CTE-I Inducible cytosolic acyl-coenzyme A thioester hydrolase LACH2 Short name= ACH2 Long chain
<b>Expression Region</b>	2-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.