



Recombinant FBJ murine osteosarcoma virus p55-v-Fos-transforming protein (V-FOS)

Product Code	CSB-EP360489FAC-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01102
Product Type	Recombinant Protein
Immunogen Species	FBJ murine osteosarcoma virus (FBJ-MSV) (Finkel-Biskis-Jenkins murine osteosarcoma virus)
Purity	≥85% (SDS-PAGE)
Sequence	MMFSGFNADY EASSFRCSSA SPAGDLSLYY HSPADSFSSM GSPVNTQDFC ADLSVSSANF IPTVTATSTS PDLQWLQPT LVSSVAPSQT RAPHYGLPT QSAGAYARAE MVKTVSGGRA QSIGRRGKVE QLSPEEEER RIRRERNKMA AAKCRNRRRE LTDTLQAETD QLEDKKSALQ TEIANLLKEK EKLEFILA AH RPACKIPDDL GFPEEMSVAS LDLTGGLPEA STPESEEAFT LPLLNDPEPK PSLEPVKSIS NVELKAEPFD DFLFPASSRP SGSETSRVSV NVDLSGSFYA ADWEPLHSNS LGMGPMVTEL EPLCTPVVTC TPLLRLPELT HAAGPVSSQR RQGSRHPDVP LPELVHYREE KHVFPQRFPS T
Source	E.coli
Target Names	V-FOS
Protein Names	Recommended name: p55-v-Fos-transforming protein
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.