



Recombinant Escherichia coli Type II methyltransferase M.EcoRI (ecoRIM)

Product Code	CSB-EP360543ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00472
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	ARNATNKLL HKAKKSKSDE FYTQYCDIEN ELQYYREHFS DKVVYCNCD PRVSNFFKYF AVNFDNLGLK KLIASCYVEN KEGFSSSEAA KNGFYEYHK ENGGKLVFDD ISVSSFCGDG DFRSSESIDL LKKS DIVTN PPFSLFREYL DQLIKYDKKF LIANVNSIT YKEVFNLIKE NKIWLGVHLG RGVSGFIVPE HYELYGTEAR IDSNGNRIIS PNNCLWLTLN DVFIRHKDLP LTRKYFGNES SYPKYDNYDA INVNKT KDIP LDYNGVMGVP ITFLHKFNPE QFELIKFRKG VDEKDLSING KCPYFRILIK NKRLQK
Source	E.coli
Target Names	ecoRIM
Protein Names	Recommended name: Modification methylase EcoRI Short name= M.EcoRI EC= 2.1.1.72 Alternative name(s): Adenine-specific methyltransferase EcoRI
Expression Region	2-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.