



Recombinant Human Haptoglobin-related protein (HPR)

Product Code	CSB-EP360576HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00739
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	L YSGNDVTDIS DDRFPKPPEI ANGYVEHLFR YQCKNYYRLR TEGDGVYTLN DKKQWINKAV GDKLPECEAV CGKPKNPANP VQRILGGHLD AKGSFPWQAK MVSHHNLTTG ATLINEQWLL TTAKNLFLNH SENATAKDIA PTLTLYVGKK QLVEIEKVVL HPNYHQVDIG LIKLKQKVLV NERVMPICLP SKNYAEVGRV GYVSGWGQSD NFKLTDHLKY VMLPVADQYD CITHYEGSTC PKWKAPKSPV GVQPILNEHT FCVGMSKYQE DTCYGDAGSA FAVHDLLEEDT WYAAGILSFD KSCAVAEYGV YVKVTSIQHW VQKTIAEN
Source	E.coli
Target Names	HPR
Protein Names	Recommended name: Haptoglobin-related protein
Expression Region	20-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.