



Recombinant Pig Pepsin A (PGA)

Product Code	CSB-YP360581PI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00791
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	≥85% (SDS-PAGE)
Sequence	I GDEPLENYLD TEYFGTIGIG TPAQDFTVIF DTGSSNLWVP SVYCSSLACS DHNQFNPDSS STFEATSQEL SITYGTGSMT GILGYDTVQV GGISDTNQIF GLSETEPGSF LYYAPFDGIL GLAYPSISAS GATPVFDNLW DQGLVSQDLF SVYLSSNDDS GSVVLLGGID SSYYTGSLNW VPVSVEGYWQ ITLDSITMDG ETIACSGGCQ AIVDTGTSLT TGPTSAIANI QSDIGASENS DGEMVISCSS IDSLPDIVFT INGVQYPLSP SAYILQDDDS CTSGFEGMDV PTSSGELWIL GDVFIRQYYT VFDRANKVG LAPVA
Source	Yeast
Target Names	PGA
Protein Names	Recommended name: Pepsin A EC= 3.4.23.1
Expression Region	60-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.