



Recombinant Mouse Renin-2 (Ren2)

Product Code	CSB-EP360582MO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00796
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SSLTDLI SPVVLNTNYLN SQYYGEIGIG TPPQTFKVIF DTGSANLWVP STKCSRLYLA CGIHSLYESS DSSSYMENGD DFTIHYGSGR VKGFLSQDSV TVGGITVTQT FGEVTELPLI PFMLAQFDGV LGMGFPAQAV GGVTPVFDHI LSQGVLKEKV FSVYYNRGPH LLGGEVVLGG SDPEHYQGDF HYVSLSKTDS WQITMKGVSV GSSTLLCEEG CEVVVDTGSS FISAPSSLK LIMQALGAKE KRLHEYVVSC SQVPTLPDIS FNLGGRAYTL SSTDYVLQYP NRRDKLCTVA LHAMDIPPT GPVWVLGATF IRKFYTEFDR HNNRIGFALA R
Source	E.coli
Target Names	Ren2
Protein Names	Recommended name: Renin-2 EC= 3.4.23.15 Alternative name(s): Angiotensinogenase Submandibular gland renin Cleaved into the following 2 chains: 1. Renin-2 heavy chain 2. Renin-2 light chain
Expression Region	64-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.