



Recombinant Woolly monkey sarcoma virus PDGF-related-transforming protein sis (V-SIS)

Product Code	CSB-MP360627WBA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01128
Product Type	Recombinant Protein
Immunogen Species	Woolly monkey sarcoma virus (WMSV) (Smian sarcoma-associated virus)
Purity	>85% (SDS-PAGE)
Sequence	MTLTWQGDPI PEELYKMLSG HSIRSFDDLQ RLLQGDSGKE DGAELDLNMT RSHSGGELES LARGKRSLSG LSVAEPAMIA ECKTRTEVFE ISRRLIDRTN ANFLVWPPCV EVQRCSGCCN NRNVQCRPTQ VQLRPVQVRK IEIVRKKPIF KKATVTLEDH LACKCEIVAA ARAVTRSPGT SQEQRAKTTQ SRVTIRTVRV RRPPKGKHRK CKHTHDKTAL KETLGA
Source	Mammalian cell
Target Names	V-SIS
Protein Names	Recommended name: PDGF-related-transforming protein sis Alternative name(s): p28sis
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.