



Recombinant Human CTD small phosphatase-like protein 2 (CTDSPL2)

Product Code	CSB-MP006145HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05D32
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRLRTRKASQ QSNQIQTQRT ARAKRKYSEV DDSLPSGG EK PSKNETGLLS SIKKFIKGST PKEERENPSK RSRIERDIDN NLITSTPRAG EKPNKQISRV RRKSQVNGEA GSYEMTNQHV KQNGKLEDNP SSGSPRRTTL LGTIFSPVFN FFSPANKNGT SGSDSPGQAV EAEIVKQLD MEQVDEITTS TTTSTNGAAY SNQAVQVRPS LNNGLEEAEE TVNRDIPPLT APVTPDSGYS SAHAEATYEE DWEVFDPPYYF IKHVPPLTEE QLNRKPALPL KTRSTPEFSL VLDLDETLVH CSLNELEDAA LTFPVLFDV IYQVYVRLRP FFREFLERMS QMYEILFTA SKKVYADKLL NILDPKKQLV RHRLFREHCV CVQGN YIKDL NILGRDLSKT IIDNSPQAF AYQLSNGIPI ESWFMDKNDN ELLKLIPFLE KLVELNEDVR PHIRDRFRLH DLLPPD
Source	Mammalian cell
Target Names	CTDSPL2
Protein Names	Recommended name: CTD small phosphatase-like protein 2 Short name= CTDSP-like 2 EC= 3.1.3.-
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.