



# Recombinant Mouse Beclin-1 (Becn1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP002658MO  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | O88597  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MEGSKASSST MQVSFVCQRC SQPLKLDTSF KILDRVTIQE LTAPLLTAAQ<br>AKPGETQEEE ANSGEFPFIE TRQDGVSRRF IPPARMMSTE SANSFTLIGE<br>ASDGGTMENL SRRLKVTGDL FDIMSGQTDV DHPLCEECDT TLLDQLDTQL<br>NVTENECQNY KRCLEILEQM NEDDSEQLQR ELKELALEEE RLIQELEDVE<br>KNRKVVAENL EKVQAEAERL DQEEAQYQRE YSEFKRQQLE LDDELKSVEN<br>QVRYAQIQLD KLKKTNVFNA TFHIWHSGQF GTINNFRLGR LPSVPVEWNE<br>INAAWGQTVL LLHALANKMG LKFQRYRLVP YGNHSYLESL TDKSKELPLY<br>CSGGLRFFWD NKFDHAMVAF LDCVQQFKEE VEKGETRFCL<br>PYRMDVEKGK IEDTGGSGGS YSIKTQFNSE EQWTKALKFM LTNLKWGLAW<br>VSSQFYNK |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | Becn1   |
| <b>Protein Names</b>     | Recommended name: Beclin-1 Alternative name(s): Coiled-coil myosin-like<br>BCL2-interacting protein   |
| <b>Expression Region</b> | 1-448   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at<br>4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the<br>contents to the bottom. Please reconstitute protein in deionized sterile water to a<br>concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final<br>concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final<br>concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,<br>storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life<br>of lyophilized form is 12 months at -20°C/-80°C.  |