



Recombinant Zea mays Actin-1 (ACT1)

Product Code	CSB-YP360866ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P02582
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	≥85% (SDS-PAGE)
Sequence	MADEDIQPIV CDNGTGMVKA GFAGDDAPRA VFPSIVGRPR HTGVMVGMGQ KDAYVGDEAQ AKRGILTLKY PIEHGIVNNW DDMENWHHTF YNELRVSPED HPVLLTEAPL NPKANREKMT QIMFETFECF AMYVAIEAVL SLYASGRRTG IVMDSGDGVS HTVPIYEGYT LPHAILRLDL AGRDLTDHLM KILTERGYSL T TSAEREIVR DIKEKLAYVA LDYEQELETA KSSSSVEKSY EMPDGQVITI GSEFRFCPEV LFQPSLVGME SPSVHEATYN SIMKCDVDIR KDLYGNVVL GGFTMFPGIA DRMSKEITSL VPSSMKVKVV APPRRKYSVW IGGSILASLS TFQQMWISKG EYDETGPGIV HMKCF
Source	Yeast
Target Names	ACT1
Protein Names	Recommended name: Actin-1
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.