



Recombinant *Saccharomyces cerevisiae* Elongation factor 1-alpha (TEF1)

Product Code	CSB-BP360931SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P02994
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGKEKSHINV VVIGHVDSGK STTTGHLIYK CGGIDKRTIE KFEKEAAELG KGSFKYAWVL DKLKAERERG ITIDIALWKF ETPKYQVTVI DAPGHRDFIK NMITGTSQAD CAILIAGGV GEFEAGISKD GQTRHALLA FTLGVRQLIV AVNKMDSVKW DESRFQEIVK ETSNFIKKVG YNPKTVPFVP ISGWNGDNMI EATTNAPWYK GWEKETKAGV VKGKTLLEAI DAIEQPSRPT DKPLRLPLQD VYKIGGIGTV PVGRVETGVI KPGMVVTFAP AGVTTEVKSVM EMHHEQLEQG VPGDNVGFNV KNVSVKEIRR GNVCGDAKND PPKGCASFNA TVIVLNHPGQ ISAGYSPVLD CHTAHIACRF DELLEKNDRR SGKKLEDHPK FLKSGDAALV KFVPSKPMCV EAFSEYPPLG RFAVRDMRQT VAVGVIKSVD KTEKAAKVTK AAQKAAKK
Source	Baculovirus
Target Names	TEF1
Protein Names	Recommended name: Elongation factor 1-alpha Short name= EF-1-alpha Alternative name(s): Eukaryotic elongation factor 1A Short name= eEF1A Translation elongation factor 1A
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.