



Recombinant Chicken Keratin, type I cytoskeletal 19 (KRT19)

Product Code	CSB-MP012534CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O93256
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MATYSFRQTT SSVAGGPCGR SLRLGGGSFR APSIHGGSGG RGVSVSSARF VSSGLGSGLG GGYGGAFSSS FSAGFGGGYG GGLGSGDGLL SGNEKTTMQN LNDRLASYLD KVRALEEANS DLETKIREWY LKQGP GPARD YSPYYKAIED LRDQILAATI DNSKVVLQID NARLAADDFK TKFETEALR MSVEADINGL RRVDELTLA RTDLELQIEN LKEELAYLKK NHEEEMSALG GQVASQVSVE VDSAPGIDLS KILADM RDQY EHMAEKNRKD AEAWFH SKTE ELNRELAVNT EQLQSSKSEV TDLRRTLQGL EIELQSQLSM KGALESTLAD TEGRYGAQLA QIQDMIGSIE AQLAELRADM ERQNSEYKML MDIKTRLEQE IATYRQLLEG QESQLFGSLS GSPDKRDKPA DGK
Source	Mammalian cell
Target Names	KRT19
Protein Names	Recommended name: Keratin, type I cytoskeletal 19 Alternative name(s): Cytokeratin-19 Short name= CK-19 GK-19 Keratin-19 Short name= K19
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Unlike its related family members, this smallest known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically expressed in the periderm, the transiently superficial layer that envelopes the developing epidermis. The type I cytokeratins are clustered in a region of chromosome 17q12-q21.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.