



# Recombinant Vesicular stomatitis Indiana virus Nucleoprotein (N)

<b>Product Code</b>	CSB-EP361079VBJ-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03521
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Vesicular stomatitis Indiana virus (strain San Juan) (VSIV)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSVTVKRIID NTVIVPKLPA NEDPVEYPAD YFRKSKEIPL YINTTKSLSD LRGYVYQGLK SGNVSIHVN SYLYGALKDI RGKLDKDWSS FGINIGKAGD TIGIFDLVSL KALDGVLPDG VSDASRTSAD DKWLPLYLLG LYRVGRTQMP EYRKKLMDGL TNQCKMINEQ FEPLVPEGRD IFDVWGND SN YTKIVA AVDM FFHMFKKHEC ASFRYGTIVS RFKDCAALAT FGHLCKITGM STEDVTTWIL NREVADEM VQ MMLPGQEIDK ADSYMPYLID FGLSSKSPYS SVKNPAFHFW GQLTALLRS TRARNARQPD DIEYTSLLTA GLLYAYAVGS SADLAQQFCV GDNKYTPDDS TGGLTTNAPP QGRDVVEWLG WFDQNRKPT PDMMQYAKRA VMSLQGLREK TIGKYAKSEF DK
<b>Source</b>	E.coli
<b>Target Names</b>	N
<b>Protein Names</b>	Recommended name: Nucleoprotein Short name= NP Alternative name(s): Nucleocapsid protein Short name= Protein N
<b>Expression Region</b>	1-422
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.