



Recombinant Rotavirus A Intermediate capsid protein VP6

Product Code	CSB-MP361082RHA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03531
Product Type	Recombinant Protein
Immunogen Species	Rotavirus A (strain RVA/SA11-Both/G3P5B[2]) (RV-A) (Simian Agent 11 (strain Both))
Purity	>85% (SDS-PAGE)
Sequence	MDVLYSLSKT LKDARDKIVE GTLYSNVSDL IQQFNQMIIT MNGNEFQTGG IGNLPIRNWN FNFGLLGTTL LNLDANYVET ARNTIDYFVD FVDNVCMDDEM VRESQRNGIA PQSDSLRKLS AIKFKRINF D NSSEYIENWN LQNRRTGF TFHKPNIFPY SASFTLNRSQ PAHDNLMGTM WLNAGSEIQV AGFDYSCAIN APANIQQFEH IVPLRRVLTT ATITLLPDAE RFSFPRVINS ADGATTWFFN PVILRPNNVE VEFLLNGQII NTYQARFGTI VARNFDTIRL SFQLMRPPNM TPAVAVLFPN AQPFEHHATV GLTLRIESAV CESVLADASE TLLANVTSVR QEYAIPVGPV FPPGMNWTDL ITNYSPSRED NLQRVFTVAS IRSMLIK
Source	Mammalian cell
Protein Names	Recommended name: Intermediate capsid protein VP6
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.