



# Recombinant Enterobacteria phage MS2 Assembly protein (A)

<b>Product Code</b>	CSB-YP361106ECZ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03610
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia phage MS2 (Bacteriophage MS2)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRAFSTLDRE NETFVPSVRV YADGETEDNS FSLKYRSNWT PGRFNSTGAK TKQWHYPSPY SRGALSVTSI DQGAYKRSGS SWGRPYEKA GFGFSLDARS CYSLFPVSQN LTYIEVPQNV ANRASTEVLQ KVTQGNFNLG VALAEARSTA SQLATQTIAL VKAYTAARRG NWRQALRYLA LNEDRKFRSK HVAGRWLELQ FGWLPLMSDI QGAYEMLT KV HLQEFLPMRA VRQVGTNIKL DGRLSYPAAN FQTTCNISRR IVIWFYINDA RLAWLSSSLGI LNPLGIVWEK VPFSFVVDWL LPVGNMLEGL TAPVGCSYMS GTVTDVITGE SIISVDAPYG WTVERQGTAK AQISAMHRGV QSVWPTTGAY VKSPFSMVHT LDALALIRQR LSR
<b>Source</b>	Yeast
<b>Target Names</b>	A
<b>Protein Names</b>	Recommended name: Assembly protein Alternative name(s): Maturation protein Short name= A protein
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.