



Recombinant Enterobacteria phage T7 Capsid assembly protein (9)

Product Code	CSB-EP361133EEB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03716
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T7 (Bacteriophage T7)
Purity	>85% (SDS-PAGE)
Sequence	MAESNADVYA SFGVNSAVMS GGSVEEHEQN MLALDVAARD GDDAIELASD EVETERDLYD NSDPFGQEDD EGRIQVRIGD GSEPTDVDTG EEGVEGTEGS EEFTPLGETP EELVAASEQL GEHEEGFQEM INIAAERGMS VETIEAIQRE YEENEELSAE SYAKLAEIGY TKAFIDSYIR GQEALVEQYV NSVIEYAGGR ERFDALYNHL ETHNPEAAQS LDNALTNRDL ATVKAINLA GESRAKAFGR KPTRSVTNRA IPAKPQATKR EGFADRSEMI KAMSDPRYRT DANYRRQVEQ KVIDSNF
Source	E.coli
Target Names	9
Protein Names	Recommended name: Capsid assembly protein Alternative name(s): Scaffolding protein
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.