



Recombinant Enterobacteria phage lambda Probable endopeptidase (K)

Product Code	CSB-YP361136ECW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03729
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Purity	≥85% (SDS-PAGE)
Sequence	MSPEDWLQAE MQGEIVALVH SHPGGLPWLS EADRRRLQVQS DLPWWLVCRG TIHKFRCVPH LTGRRFEHGV TDCYTLFRDA YHLAGIEMPD FHREDDWWRN GQNLVLDNLE ATGLYQVPLS AAQPGDVLLC CFGSSVPNHA AIYCGDGELL HHIPEQLSKR ERYTDKWQRR THSLWRHRAW RASAFTGIYN DLVAASTFV
Source	Yeast
Target Names	K
Protein Names	Recommended name: Probable endopeptidase EC= 3.4.-.- Alternative name(s): Tail assembly protein K
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.