



Recombinant Enterobacteria phage lambda Uncharacterized nin region protein ORF290

Product Code	CSB-MP361147ECW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03766
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Purity	≥85% (SDS-PAGE)
Sequence	MINVVSFSGG RTSAYLLWLM EQKRRAGKDV HYVFM DTGCE HPMTYRFVRE VVKFWDIPLT VLQVDINPEL GQPNGYTVWE PKDIQTRMPV LKPFIDMVKK YGTPYVGGAF CTDRLKLVPF TKYCDDHFGR GNYTTWIGIR ADEPKRLKPK PGIRYLAELS DFEKEDILAW WKQQPFDLQI PEHLGNCIFC IKKSTQKIGL ACKDEEGLQR VFNEVITGSH VRDGHRETPK EIMYRGRMSL DGI AKMYSEN DYQALYQDMV RAKRFDTGSC SESCEIFGGQ LDFDFGREAA
Source	Mammalian cell
Protein Names	Recommended name: Uncharacterized nin region protein ORF290
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.