



Recombinant Enterobacteria phage lambda Serine/threonine-protein phosphatase

Product Code	CSB-MP361149ECW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03772
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Purity	>85% (SDS-PAGE)
Sequence	MRYYEKIDGS KYRNIWVVGD LHGCYTNLMN KLDTIGFDNK KDLLISVGD VDRGAENVEC LELITFPWFR AVRGNHEQMM IDGLSERGNV NHWLLNGGGW FFNLDYDKEI LAKALAHKAD ELPLIIELVS KDKKYVICHA DYPFDEYEFK KPVDPHQVIV NRERISNSQN GIVKEIKGAD TFIFGHTPAV KPLKFANQMY IDTGAVFCGN LTLIQVQGEG A
Source	Mammalian cell
Protein Names	Recommended name: Serine/threonine-protein phosphatase EC= 3.1.3.16
Expression Region	1-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.