



Recombinant *Saccharomyces cerevisiae* Intron-encoded RNA maturase bl4 (bl4), partial

| | |
|--------------------------|---|
| Product Code | CSB-EP361174SVG |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P03879 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Purity | ≥85% (SDS-PAGE) |
| Source | <i>E.coli</i> |
| Target Names | Bl4 |
| Protein Names | Recommended name: Intron-encoded RNA maturase bl4 Alternative name(s): RNA maturase SCbl4 Cleaved into the following 2 chains: 1. Truncated, nonfunctional cytochrome b 2. RNA maturase bl4 EC= 3. 3.1.-.- |
| Expression Region | - |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |