



Recombinant Human papillomavirus type 11 Regulatory protein E2 (E2)

Product Code	CSB-YP361192HMG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04015
Product Type	Recombinant Protein
Immunogen Species	Human papillomavirus type 11
Purity	>85% (SDS-PAGE)
Sequence	MEAIKRLDA CQDQLELYE ENSIDIHKHI MHWKCIRLES VLLHKAKQMG LSHIGLQVVP PLTVSETKGH NAIEMQMHLE SLAKTQYGVE PWTLQDTSYE MWLTPPKRCF KKQGNTVEVK FDGCEDNVME YVWVTHIYLQ DNDSWVKVTS SVDAGIYYT CGQFKTYVFN FNKEAQKYGS TNHWEVCYGS TVICSPASVS STVREVSIAE PTTYTPAQT APTVSACTTE DGVSAPPRKR ARGPSTNNTL CVANIRSVDS TINNIVTDNY NKHQRRNNCH SAATPIVQLQ GDSNCLKCFR YRLNDKYKHL FELASSTWHW ASPEAPHKNA IVTLTYSSEE QRQQFLNSVK IPPTIRHKVG FMSLHLL
Source	Yeast
Target Names	E2
Protein Names	Recommended name: Regulatory protein E2
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.