



Recombinant Human Nucleolar protein 4 (NOL4)

Product Code	CSB-MP015922HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94818
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESERDMYRQ FQDWCLRTYG DSGKTKTVTR KKYERIVQLL NGSESSSTDN AKFKFWVSKS GFQLGQPDEV RGGGGGAKQV LYVPVKTTDG VGVDEKLSLR RVAVVEDFFD IIYSMHVETG PNGEQIRKHA GQKRTYKAIS ESYAFLPREA VTRFLMSCSE CQKRMHLNPD GTDHKDNGKP PTLVTSMIDY NMPITMAYMK HMKLQLLSQ QDEDESSIES DEFDMSDSTR MSAVNSDLSS NLEERMQSPQ NLHGQQDDDS AAESFNGNET LGHSSIASGG THSREMGSND SDGKTGLEQD EQPLNLSDSP LSAQLTSEYR IDDHNSNGKN KYKNLLISDL KMEREARENG SKSPAHSYSS YDSGKNESVD RGAEDLSLNR GDEDEDDHED HDDSEKVNET DGVEAERLKA FNMFVRLFVD ENLDRMVPIS KQPKEKIQAI IDSCRRQFPE YQERARKRIR TYLKSCRRMK RSGFEMSRPI PSHLTSAVAE SILASACESE SRNAAKRMRL ERQQDESAPA DKQCKPEATQ ATYSTSAVPG SQDVLYINGN GTYSYHSYRG LGGGLLNLD ASSSGPTDLS MKRQLATSSG SSSSSNSRPQ LSPTEINAVR QLVAGYRESA AFLRSADEL ENLILQQN
Source	Mammalian cell
Target Names	NOL4
Protein Names	Recommended name: Nucleolar protein 4 Alternative name(s): Nucleolar-localized protein
Expression Region	1-638
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.