



Recombinant Human F-BAR and double SH3 domains protein 2 (FCHSD2)

Product Code	CSB-EP008549HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O94868
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MQPPPRKVKV TQELKNIQVE QMTKLQAKHQ AECDLLEDMR TFSQKKAIE REYAQGMQKL ASQYLKRDWP GVKADDRNDY RSMYPVWKSF LEGTMQVAQS RMNICENYKN FISEPARTVR SLKEQQLKRC VDQLTKIQTE LQETVKDLAK GKKKYFETE MAHAVREKAD IEAKSKLSLF QSRISLQKAS VKLKARRSEC NSKATHARND YLLTAAANA HQDRYYQTDL VNIMKALDGN VYDHLKDYLI AFSRTELETC QAVQNTFQFL LENS SKVVRD YNLQLFLQEN AVFHKPQPFQ FQPCSDTSR QLESETGTTE EHSLNKEARK WATRVAREHK NIVHQQRVLN DLECHGAAVS EQSRAELEQK IDEARENIRK AEI IKLKAEA RLDLLKQIGV SVDTWLKSAM NQVMEELENE RWARPPAVTS NGTLHSLNAD TERE EEEFE DNMDVFDDSS SPSGTLRNY PLTCKVVYSY KASQPDELTI EEHEVLEVIE DGMEDWVKA RNKVGQVGYV PEKYLQFPTS NSLLSMLQSL AALDSRSHS SNSTEAE LVS GSLNGDASVC FVKALYDYEG QTDELSPFE GAIIRILNKE NQDDDFWEG EFNGRIGVFP SVLVEELSAS ENGDTPWMRE IQISPSPKPH ASLPPLPLYD QPPSSPYSP DKRSSLYFPR SPSANEKSLH AESPGFSQAS RHTPETS YGK LRPVRAAPP PTQNHRRPAE KIEDVEITLV </p>
Source	E.coli
Target Names	FCHSD2
Protein Names	Recommended name: FCH and double SH3 domains protein 2 Alternative name(s): Carom SH3 multiple domains protein 3
Expression Region	1-740
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.