



# Recombinant Human herpesvirus 1 Transcriptional regulator ICP22 (US1)

<b>Product Code</b>	CSB-BP361276HWY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P04485
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADISPGAFA PCVKARRPAL RSPPLGTRKR KRPSRPLSSE SEVESDTALE SEVESETASD STESGDQDEA PRIGGRRAPR RLGGRFFLDM SAESTTGTET DASVSDDPDD TSDWSYDDIP PRPKRARVNL RLTSPPDRRD GVIFPKMGRV RSTRETQPRA PTPSAPSPNA MLRRSVRQAQ RRSSARWTPD LGYMRQCINQ LFRVLRVARD PHGSANLRH LIRDCYLMGY CRARLAPRTW CRLQLVSGGT WGMHLRNTIR EVEARFDATA EPVCKLPCLE TRRYGPECDL SNLEIHLSAT SDDEISDATD LEAAGSDHTL ASQSDTEDAP SPVTLETPEP RGLAVRLED EFGFDWTPQ EGSQPWLSAV VADTSSVERP GPSDSGAGRA AEDRKCLDGC RKMRFSTACP YPCSDTFLRP
<b>Source</b>	Baculovirus
<b>Target Names</b>	ICP22
<b>Protein Names</b>	Recommended name: Transcriptional regulator ICP22 Alternative name(s): Immediate-early protein IE68 Infected cell protein 22 Short name= ICP22
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.