




# Recombinant Mouse Receptor-type tyrosine-protein phosphatase N2 (Ptpn2), partial

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP019060MO                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Relevance</b>         | Implicated in development of nervous syst and pancreatic endocrine cells.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Abbreviation</b>      | Recombinant Mouse Ptpn2 protein, partial                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                  |
| <b>Uniprot No.</b>       | P80560                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Alias</b>             | PTP IA-2betaProtein tyrosine phosphatase-NP ;PTP-NP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Purity</b>            | ≥ 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Sequence</b>          | RGRQLPGRLGCLFEDGLCGSLETCVNDGVFGRCQKVPVMDTYRYEVPPGAL<br>LHLKVTLQKLSRTGFTWQDDYTQRVIAQELANLPKAYLWHGEASGPAPSLQQ<br>NADNEKWFSLEREVALAKTLRRYLPYLELLSQPTANAHSRIDHETRPAGKED<br>SSPENILTYVAHTSALTYPATRAKYPDNLLRPFSRLQPDELSPKVDGDIDKQK<br>LIAALGAYTAQRLPGENDPEPRYL VHGSARAPRPFSATALSQRWPPPPGDAK<br>DSPSMDDDTLLQSLKDLQNSEVDRLGPLKEEKADSVAGAIQSDPAEGSQE<br>SHGRGAEGQPREQTDAPETMLQDHRLSEVDDPVYKEVNRLSFQLGDLKDY<br>GSPLLPEGPLLEKSSREEIKKSEQPEEVLSSSEETAGVEHVRSTYSKDLFERK<br>PNSEPQPRRLEDQFQNRAPELWEDEESLKLAQGGPSGGLQLEVQPSEEQQ<br>GYILTGNNPLSPEKGKQLMDQVAHILRVPSSFFADIKVLGPAVTFKVSANIQNM<br>TTADVIKAAADNKDQLEKATGLTILQSGIRPKGKLLKLLPHQEEQEDSTK |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Target Names</b>      | Ptpn2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | PTP IA-2betaProtein tyrosine phosphatase-NP ;PTP-NP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Expression Region</b> | 28-600aa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tag Info</b>          | N-terminal 6xHis-tagged                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Mol. Weight</b>       | 65.7 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | Extracellular Domain                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.