



# Recombinant Human RING finger protein 37 (UBOX5)

<b>Product Code</b>	CSB-MP025507HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O94941
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVINLCLPQF RPRIHCNKIS ADGYEVENLI SEDLTKRSHG FRTEYFIKPP VYVTVSFPFN VEICRINIDL TAGGGQNVTV LEMYSASSS RVSWNTPQCR TLGPAEPSVP DKEAFTLVGK VLLKNQSQVV FSHRGFKARP PFGAMEATLP SPAVVAQELW NKGALSLSHV AHLRICITHV TGGGIPCIKR LEVWGQPAKT CSQEVIDSIL LVTSENLPQD VALQAPALPM ESDCDPGDQP ESQQAPSSLQ KLAEIIQDVP EEFLDPITLE IMPCPMLLPS GKVIDQSTLE KCNRSEATWG RVPSDPFTGV AFTPHSQPLP HPSLKARIDH FLLQHSIPGC HLLGRAQTAL AVIPSSIVLP SQKRKIEQAE HVPDSNFGVN ASCFSATSPL VLPTTSEHTA KKMKATNEPS LTHMDCSTGP LSHEQKLSQS LEIALASTLG SMPSTARLT RGQLQHLGTR GSNTSWRPGT GSEQPGSILG PECASCKRVF SPYFKKEPVY QLPCGHLLCR PCLGEKQRSL PMTCTACQRP VASQDVLRVH F
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UBOX5
<b>Protein Names</b>	Recommended name: RING finger protein 37 Alternative name(s): U-box domain-containing protein 5 Ubiquitin-conjugating enzyme 7-interacting protein 5
<b>Expression Region</b>	1-541
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.