



Recombinant Human RING finger protein 37 (UBOX5)

Product Code	CSB-YP025507HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94941
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVINLCLPQF RPRIHCNKIS ADGYEVENLI SEDLTKRSHG FRTEYFIKPP VYVTVSFPFN VEICRINIDL TAGGGQNVTV LEMYSASSS RVSWNTPQCR TLGPAEPSVP DKEAFTLVGK VLLKNQSQVV FSHRGFKARP PFGAMEATLP SPAVVAQELW NKGALSLSHV AHLRICITHV TGGGIPCIKR LEVWGQPAKT CSQEVIDSIL LVTSENLPQD VALQAPALPM ESDCDPGDQP ESQQAPSSLQ KLAEIIQDVP EEFLDPITLE IMPCPMLLPS GKVIDQSTLE KCNRSEATWG RVSPDPFTGV AFTPHSQPLP HPSLKARIDH FLLQHSIPGC HLLGRAQTAL AVIPSSIVLP SQKRKIEQAE HVPDSNFGVN ASCFSATSPL VLPTTSEHTA KKMKATNEPS LTHMDCSTGP LSHEQKLSQS LEIALASTLG SMPSTARLT RGQLQHLGTR GSNTSWRPGT GSEQPGSILG PECASCKRVF SPYFKKEPVY QLPCGHLLCR PCLGEKQRSL PMTCTACQRP VASQDVLRVH F
Source	Yeast
Target Names	UBOX5
Protein Names	Recommended name: RING finger protein 37 Alternative name(s): U-box domain-containing protein 5 Ubiquitin-conjugating enzyme 7-interacting protein 5
Expression Region	1-541
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.