



Recombinant Pseudomonas aeruginosa Type IV major pilin protein PilA (pilA), partial

Product Code	CSB-MP361358EZC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04739
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / PAO1 / 1C / PRS 101 / LMG 12228)
Purity	>85% (SDS-PAGE)
Sequence	PQYQNYVARSEGASALATINPLKTTVEESLSRGIAGSKIKIGTTASTATETYVGV EPDANKLGVIAVAIEDSGAGDITFTFQTGTSSPKNATKVITLNRTADGVWACKS TQDPMFTPKGCDN
Source	Mammalian cell
Target Names	pilA
Protein Names	Recommended name: Fimbrial protein Alternative name(s): Pilin Strain PAO
Expression Region	28-149aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.