




# Recombinant Mouse Serum amyloid A-3 protein (Saa3), partial

<b>Product Code</b>	CSB-YP361411MO
<b>Relevance</b>	Major acute phase reactant. Apolipoprotein of the HDL complex.
<b>Abbreviation</b>	Recombinant Mouse Saa3 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P04918
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	RWVQFMKEAGQGSRDMWRAYSDMKKANWKNSDKYFHARGNYDAARRGPG GAWAAKVISDAREAVQKFTGHGAEDSRADQFANEWGRSGKDPNHFRPAGLP KRY
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	Saa3
<b>Protein Names</b>	Recommended name: Serum amyloid A-3 protein
<b>Expression Region</b>	20-122aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	13.8kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.