



Recombinant Bordetella pertussis Pertussis toxin subunit 3 (ptxC)

Product Code	CSB-MP361434BUA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04979
Product Type	Recombinant Protein
Immunogen Species	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
Purity	≥85% (SDS-PAGE)
Sequence	VA PGIVIPPKAL FTQQGGAYGR CPNGTRALTV AELRGNAELQ TYLRQITPGW SIYGLYDGT Y LGQAYGGIIK DAPPGAGFIY RETFCITTIY KTGQPAADHY YSKVTATRL ASTNSRLCAV FVRDGQSVIG ACASPYEGRY RDMYDALRRL LYMIYMSG LA VRVHVSKEEQ YYDYEDATFQ TYALTGISLC NPAASIC
Source	Mammalian cell
Target Names	ptxC
Protein Names	Recommended name: Pertussis toxin subunit 3 Short name= PTX S3 Alternative name(s): Islet-activating protein S3 Short name= IAP S3
Expression Region	29-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.