



Recombinant Hantaan virus Nucleoprotein (N)

Product Code	CSB-MP361470HCB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05133
Product Type	Recombinant Protein
Immunogen Species	Hantaan virus (strain 76-118) (Korean hemorrhagic fever virus)
Purity	>85% (SDS-PAGE)
Sequence	MATMEELQRE INAHEGQLVI ARQKVRDAEK QYEKDPDELN KRTLTDREGV AVSIQAKIDE LKRQLADRIA TGKNLGKEQD PTGVEPGDHL KERSMLSYGN VLDLNHLDID EPTGQTADWL SIIVYLTSFV VPILLKALYM LTTRGRQTTK DNKGTRIRFK DDSSFEDVNG IRKPKHLYVS LPNAQSSMKA EEITPGRYRT AVCGLYPAQI KARQMISPVM SVIGFLALAK DWSDRIEQWL IEPCKLLPDT AAVSLGGPA TNRDYLRQRQ VALGNMETKE SKAIRQHAEA AGCSMIEDIE SPSSIWVFAG APDRCPPTCL FIAGIAELGA FFSILQDMRN TIMASKTVGT SEEKLRKKSS FYQSYLRRTQ SMGIQLGQRI IVLFMVAWGK EAVDNFHLGD DMDPELRTLA QSLIDVKVKE ISNQEPLKL
Source	Mammalian cell
Target Names	N
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= Protein N
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.