



Recombinant Avian sarcoma virus Viral jun-transforming protein (JUN)

Product Code	CSB-MP361514ATE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05411
Product Type	Recombinant Protein
Immunogen Species	Avian sarcoma virus (strain 17)
Purity	>85% (SDS-PAGE)
Sequence	MSAKMEPTFY EDALNASFAP PESGGYGYNN ADILTSPDVG LLKLSPELE RLIIQSSNGL ITTPTPTQF LCPKNVTDEQ EGFAEGFVRA LAELHNQNTL PSVTSAAQPV SGGMAPVSSM AGGGSFNTSL HSEPPVYANL SNFNPNALNS APNYNANRMG YAPQHHINPQ MPVQHPRQLA LKEEPQTVPE MPGETPPLFP IDMESQERIK AERKRMRNRI AASKSRKRKL ERIARLEEKV KTLKAQNSEL ASTANMLREQ VAQLKQKVMN HVNSGCQLML TQQLQTF
Source	Mammalian cell
Target Names	JUN
Protein Names	Recommended name: Viral jun-transforming protein Short name= v-Jun
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.