



Recombinant Escherichia coli Glucitol/sorbitol-specific phosphotransferase enzyme IIA component (srIB)

Product Code	CSB-EP361600ENV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05706
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTVIYQTTIT RIGASAI DAL SDQMLITFRE GAPADLEEYC FIHCHGELKG ALHPGLQFSL GQHRYPVTAV GSVAEDNLRE LGHVTLRF DG LNEAEFPGTV HVAGPVPDDI APGSVLKFES VKE
Source	E.coli
Target Names	srIB
Protein Names	Recommended name: Glucitol/sorbitol-specific phosphotransferase enzyme IIA component EC= 2.7.1.- Alternative name(s): EIIA-Gut EIII-Gut PTS system glucitol/sorbitol-specific EIIA component
Expression Region	1-123
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.