



Recombinant Simian immunodeficiency virus agm.vervet Envelope glycoprotein gp160 (env), partial

Product Code	CSB-EP361658SHH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05886
Product Type	Recombinant Protein
Immunogen Species	Simian immunodeficiency virus agm.vervet (isolate AGM TYO-1) (SIV-agm.ver) (Simian immunodeficiency virus African green monkey vervet)
Purity	>85% (SDS-PAGE)
Source	E.coli
Target Names	env
Protein Names	Recommended name: Envelope glycoprotein gp160 Alternative name(s): Env polyprotein Cleaved into the following 2 chains: 1. Surface protein gp120 Short name= 2. SU Alternative name(s): Glycoprotein 120 Short name= gp120 Tr
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.