



Recombinant Simian immunodeficiency virus agm.vervet Virion infectivity factor (vif)

Product Code	CSB-MP361667SHH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05904
Product Type	Recombinant Protein
Immunogen Species	Simian immunodeficiency virus agm.vervet (isolate AGM TYO-1) (SIV-agm.ver) (Simian immunodeficiency virus African green monkey vervet)
Purity	>85% (SDS-PAGE)
Sequence	MNPNKEWVMR VTWKVPGDLI TKWQGIVRYW MRQRNLKWNV YMHYQITWAW YTMSRYVIPI GKHGEICVDL YWHLTPEQGW LSTYAVGIQY VSNLESKYRT ELDPATADSI IHGHYFNCFK ERAIQQALRG HRFVFCQFPE GHKSTGQVPS LQYLALLAHQ NGLRERSKRG KTRRSRNLGS KQGAVGQMAK RYVTRSQPGG EAAFWERTPV PSMELLSGGR RKTWYSHDGK GLQIL
Source	Mammalian cell
Target Names	vif
Protein Names	Recommended name: Virion infectivity factor Short name= Vif Alternative name(s): Q protein SOR protein
Expression Region	1-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.