



Recombinant Human Putative protein PLEKHA9 (PLEKHA8P1)

Product Code	CSB-EP018163HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95397
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSELRLCCDL LVQQVDKTK E VTTTGVSNSE EGIDVGTLK STCNTFLKTL EECMQIANAA FTSELLYHTP PGSPQLAMLK SSKMKHPIIP IHNSLERQTE LSTCENGSLN MEINGEEEIL MKNKNSLYLK SAEIDCSISS EENTDDNITV QGEIMKEDRM ENLKNHDNNL SQSGSDSSCS PECLWEEGKE VIPTFFSTMN TSFSDIELLE DSGIPTEAFL ASCCAVVPVL DKLGPVTFAP VKMDLVENIK KVNQKYITNK EEFTTLQKIV LHEVEADVAQ VRNSATEALL WLKRGLKFLK GFLTEVKNGE KDIQTALNNA YGKTLRQHHG WVVRGVFALA LRATPSYEDF VAALTVKEGD HRKEAFSIGM QRDL SLYLPA MKKQMAILDA L
Source	E.coli
Target Names	PLEKHA8P1
Protein Names	Recommended name: Putative protein PLEKHA9 Alternative name(s): Pleckstrin homology domain-containing family A member 8 pseudogene 1
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.