



Recombinant Salmonella typhimurium Ketol-acid reductoisomerase (ilvC)

Product Code	CSB-MP361693SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P05989
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	<p> ANYFNTLNL RQQLAQLGKC RFMGRDEFAD GASYLQGKKV VIVGCGAQGL NQGLNMRDSG LDISYALRKE AIAEKRASWR KATENGFKVG TYEELIPQAD LVVNLTPDKQ HSDVVRVSVQP LMKDGAALGY SHGFNIVEVG EQIRKDITVV MVAPKCPGTE VREEYKRGFG VPTLIAVHPE NDPQGEGMAI AKAWAAATGG HRAGVLESSF VAEVKSDLMG EQTILCGMLQ AGSLLCFDKL VAEGTDPAYA EKLIQFGWET ITEALKQGGI TLMMDRLSNP AKLRAYALSE QLKEIMAPLF QKHMDDIISG EFSSGMMADW ANDDKLLTW REETGKTAFE TAPQYEGKIG EQEYFDKGV L MIAMVKAGVE LAFETMVDSG IIEESAYYES LHELPLIANT IARKRLYEMN VVISDTAEYG NYLFSYACVP LLKPFMAELQ PGDLGSAIPE GAVDNAQLRD VNDAIRSHAI EQVGKCLRGY MTDMKRIAVA G </p>
Source	Mammalian cell
Target Names	ilvC
Protein Names	Recommended name: Ketol-acid reductoisomerase EC= 1.1.1.86 Alternative name(s): Acetohydroxy-acid isomeroreductase Alpha-keto-beta-hydroxylacil reductoisomerase
Expression Region	2-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.